



PACIFIC
BIOSCIENCES®



Using the Iso-Seq Application on SMRT Link and BioConda

Elizabeth Tseng, PacBio

 @Magdoll



Why use Iso-Seq analysis?

ISO-SEQ ANALYSIS MAIN FEATURES

- No reference genome required
- No transcriptome assembly required
- Recovers full-length (5' to 3') transcripts
- Yields highly accurate (>99%) transcripts

ISOSEQ.HOW

Iso-Seq Docs

[Iso-Seq Home](#)

[Nomenclature](#)

[Single cell](#) ▾


[Clustering](#) ▾

[General FAQ](#)

[Changelog](#)

This site uses [Just the Docs](#), a documentation theme for Jekyll.

🔍 Search Iso-Seq Docs [File an issue](#)



Iso-Seq


Scalable
De Novo
Isoform Discovery
from PacBio HiFi Reads

Iso-Seq contains the newest tools to identify transcripts in PacBio single-molecule sequencing data. Starting in SMRT Link v6.0.0, those tools power the *Iso-Seq GUI-based analysis* application. A composable workflow of existing tools and algorithms, combined with a new clustering technique, allows to process the ever-increasing yield of PacBio. Starting with version 3.4, support for UMI and cell barcode based deduplication has been added.

Availability

Latest version can be installed via bioconda package `isoseq3`.


HIFI READS FROM CCS




FULL-LENGTH READS HAVE 5' AND 3' PRIMERS



HiFi reads




FL reads



REMOVE CONCATEMERS AND POLY(A) TAILS




FL reads



FLNC reads



CLUSTER TO GET ISOFORMS



- High Quality (HQ):
accuracy $\geq 99\%$ and ≥ 2 FLNC read support
- Low Quality (LQ):
accuracy $< 99\%$ and ≥ 2 FLNC read support


MAP AND COLLAPSE ISOFORMS



High-Quality, Full-Length Polished Isoforms



Map to Reference Genome
Minimap2 (pbmm2)



BENEFITS OF ISO-SEQ ANALYSIS APPLICATION

- High-quality transcripts
- Full-Length Non-concatemer reads
- Mapped & collapsed isoforms
- Removes artifacts
- Removes poly(A) tails



Iso-Seq Analysis Using pbBioConda

INSTRUCTIONS TUTORIAL






Follow the instructions tutorial for installing all the software needed.

- If you do not have an HPC server to install pbbioconda, you should have already:
 - Create an AWS account
 - Create an AWS Linux Instance to run Iso-Seq 3 Analysis Pipeline
 - Connect to your AWS Instance
- Upgrades and Install Software

DOWNLOAD THE DATA

https://downloads.pacbcloud.com/public/dataset/ISMB_workshop/

Index of /public/dataset/ISMB_workshop/iseq3

Name	Last modified	Size	Description
 Parent Directory		-	
 results/	2020-09-23 07:31	-	
 alz.ccs.bam	2020-06-15 11:52	84M	
 iseq_primers.fasta	2020-09-23 07:03	62	
 run.sh	2020-09-23 07:23	430	

Example:

```
$ wget -nv https://downloads.pacbcloud.com/public/dataset/ISMB_workshop/iseq3/alz.ccs.bam
```

SPECIFY ISO-SEQ PRIMERS


```
$ more primers.fasta
```

```
>5p
```

```
GCAATGAAGTCGCAGGGTTGGG
```

```
>3p
```

```
GTACTCTGCGTTGATACCACTGCTT
```



INPUT CCS BAM FILE

```
$ samtools view -h alz.ccs.bam
```

```
m141008_060349_42194_c100704972550000001823137703241586_s1_p0/63/ccs4*0255  
**00CCCGGGGATCCTCTAGAATGC~~~~~RG:Z:83ba013f np:i:35  
rq:f:0.999682 sn:B:f,11.3175,6.64119,11.6261,14.5199 zm:i:63
```

REFERENCE GENOME

```
$ grep '>' hg38.fa # to list the headers per chromosome
```

```
>chr1 AC:CM000663.2 gi:568336023 LN:248956422 rl:Chromosome M5:6aef897c3d6ff0c78a
ff06ac189178dd AS:GRCh38
>chr2 AC:CM000664.2 gi:568336022 LN:242193529 rl:Chromosome M5:f98db672eb0993dcfd
abafe2a882905c AS:GRCh38
>chr3 AC:CM000665.2 gi:568336021 LN:198295559 rl:Chromosome M5:76635a41ea913a405d
ed820447d067b0 AS:GRCh38
>chr4 AC:CM000666.2 gi:568336020 LN:190214555 rl:Chromosome M5:3210fecf1eb92d5489
da4346b3fddc6e AS:GRCh38
>chr5 AC:CM000667.2 gi:568336019 LN:181538259 rl:Chromosome M5:a811b3dc9fe66af729
dc0dddf7fa4f13 AS:GRCh38 hm:47309185-49591369
...
```


SOFTWARE INSTALLATION CHECK

Access to your conda environment

```
$ source activate <name of your environment>
```

Check your installation

```
$ isoseq3 --version
```

```
isoseq3 3.4.x
```

```
$ lima --version
```

```
lima 1.11.0
```

```
$ pbmm2 --version
```

```
pbmm2 1.3.0
```




ISO-SEQ WORKFLOW

Clustering / High-level Workflow


High-level Workflow

The high-level workflow depicts files and processes:




← Most SQ2 projects will directly give you ccs.bam now


ISO-SEQ WORKFLOW




- Use **polished** CCS reads
- Only full-pass ZMWs




- Barcoded and unbarcoded cDNA primer removal
- Orientation
- Unwanted primer combination removal



- PolyA tail trimming
- Concatemer removal





- Hierarchical, $n \cdot \log(n)$ clustering, alignment of shorter to longer sequences
- Iterative cluster merging
- Generate consensus for each read cluster using QV guided PoA



- One consensus per read cluster

- Align to reference genome
- Remove redundancy





PRIMER REMOVAL & DEMULTIPLEXING

Command line:

```
lima --isoseq --dump-clips --peek-guess -j 24\
alz.ccs.bam isoseq_primers.fasta alz.demult.bam
```

Input files:

```
alz.ccs.bam #HiFi reads
isoseq_primers.fasta #Iso-Seq primers
```





Output files:

```
alz.demult.bam
```

Options:

```
--isoseq #specialized isoseq option for lima
--dump-clips # show the clipped primers
--peek-guess # remove spurious false positive signal
-j 24 # Number of threads to use
```

PRIMER REMOVAL & DEMULTIPLEXING




After completion, you will see the following files:

```
$ ls -ltrh
```

```
alz.demult.json
alz.demult.lima.clips
alz.demult.lima.counts
alz.demult.lima.guess
alz.demult.lima.report
alz.demult.lima.summary
alz.demult.5p--3p.bam
alz.demult.5p--3p.bam.pbi
alz.demult.5p--3p.subreadset.xml
```

#lima reports



TRIMMING POLY(A) TAILS AND CONCATEMER REMOVAL

Command line:

```
isoseq3 refine --require-polya\  
alz.demult.5p--3p.bam\ isoseq_primers.fasta  
alz.flnc.bam
```

Input files:

```
alz.demult.5p--3p.bam  
isoseq_primers.fasta
```





Output files:

```
alz.flnc.bam
```

Options:

```
--require-polya #if your transcripts have a polyA tail
```

TRIMMING POLY(A) TAILS AND CONCATEMER REMOVAL



After completion, you will see the following files:

```
$ ls -ltrh
```

alz.flnc.bam

alz.flnc.bam.pbi

alz.flnc.consensusreadset.xml


alz.flnc.filter_summary.json

alz.flnc.report.csv

#isoseq3 refine reports



ISOFORMS



Command line:

```
isoseq3 cluster alz.flnc.bam alz.polished.bam \
--verbose --use-qvs
```

Input files:

alz.flnc.bam


Output files:

alz.polished.bam

Options:

--verbose #if your transcripts have a polyA tail
 --use-qvs #Use CCS QVs, sets --poa-cov 100

ISOFORMS



After completion, you will see the following files:

```
$ ls -ltrh
```

alz.polished.bam

alz.polished.bam.pbi

alz.polished.transcriptset.xml

alz.polished.cluster

alz.polished.cluster_report.csv

alz.polished.hq.bam

alz.polished.hq.bam.pbi

alz.polished.hq.fasta.gz

alz.polished.lq.bam

alz.polished.lq.bam.pbi

alz.polished.lq.fasta.gz


Because the ccs input is Polished, the isoseq3 cluster output is already polished!

#isoseq3 cluster reports

#high quality isoforms (≥ 0.99)

#low quality isoforms (< 0.99)

MAP



Command line:

```
pbmm2 align hg38.fa alz.polished.hq.bam
alz.aligned.bam
-j 24 --preset ISOSEQ -sort --log-level INFO
```

Input files:

```
alz.polished.hq.bam
hg38.fa
```


Output files:

```
alz.aligned.bam
```

Options:

```
-j 24 #Number of threads to use
--preset ISOSEQ #select the alignment mode
--sort #Generate sorted BAM file
--log-level INFO #show progress
```

MAP




After completion, you will see the following files:

```
$ ls -ltrh
```

```
alz.aligned.bam
alz.aligned.bam.bai
```

COLLAPSE



Command line:

```
isoseq3 collapse alz.aligned.bam alz.collapsed.gff
```


Input files:

alz.aligned.bam

Output files:

alz.collapsed.gff

COLLAPSE



After completion, you will see the following files:

```
$ ls -ltrh
```

```
alz.collapsed.report.json
alz.collapsed.abundance.txt
alz.collapsed.read_stat.txt
alz.collapsed.group.txt
alz.collapsed.gff
alz.collapsed.fasta
```

#report, stats and list

ISO-SEQ ANALYSIS TERMINOLOGY

NAME	ABBR	EXPLANATION
Full-Length Reads	FL Reads	CCS reads with 5' and 3' cDNA primers removed
Full-Length, Non-Concatemer Reads	FLNC Reads	CCS reads with 5' and 3' cDNA primers, polyA tail, and concatemers removed
High-Quality Isoforms	HQ Isoforms	Polished transcript sequences with predicted accuracy $\geq 99\%$ & ≥ 2 FLNC
Low-Quality Isoforms	LQ Isoforms	Polished transcript sequences with predicted accuracy $< 99\%$ & ≥ 2 FLNC



PACBIO®

www.pacb.com

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. BluePippin and SageELF are trademarks of Sage Science. NGS-go and NGSengine are trademarks of GenDx. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.