

# Long-read RNA-seq transcriptomic analysis with TALON and Swan




Fairlie Reese  
UC Davis IsoSeq Workshop  
September 1, 2021




TALON identifies, quantifies, and filters full-length transcript isoforms

# TALON identifies, quantifies, and filters full-length transcript isoforms




- Preprocessing with TranscriptClean corrects long-read sequencing artifacts

# TALON identifies, quantifies, and filters full-length transcript isoforms




- Preprocessing with TranscriptClean corrects long-read sequencing artifacts
- TALON annotates reads to transcript models
- Assigns novel models an informative novelty type

# TALON identifies, quantifies, and filters full-length transcript isoforms





- Preprocessing with TranscriptClean corrects long-read sequencing artifacts
- TALON annotates reads to transcript models
- Assigns novel models an informative novelty type
- Filters novel transcript models for reproducibility and internal priming

# TALON identifies, quantifies, and filters full-length transcript isoforms




- Preprocessing with TranscriptClean corrects long-read sequencing artifacts
- TALON annotates reads to transcript models
- Assigns novel models an informative novelty type
- Filters novel transcript models for reproducibility and internal priming
- Quantifies transcript abundances

# TALON identifies, quantifies, and filters full-length transcript isoforms



# Transcript isoforms differences are difficult to visually interpret and analyze


- Difficult to tell what the difference between any two transcript models is
- Current transcriptome analysis platforms do not prioritize isoform-level differences from long-read data
- How do we compare long-read transcriptomes between samples?





# Swan is a Python library designed for the analysis and visualization of full-length transcriptomes

- Differential expression tests on the gene and transcript level
- Isoform/TSS/TES switching detection
- Novel exon skipping and intron retention event detection
- Unique and interpretable visualizations and reports with a robust set of options
- Compatibility with single-cell data




Swan represents transcript models as a series of splice sites and intronic or exonic connections between them



**TALON DB  
or GTF**


Swan represents transcript models as a series of splice sites and intronic or exonic connections between them

- Assign each unique splice site to a node



Swan represents transcript models as a series of splice sites and intronic or exonic connections between them


- Assign each unique splice site to a node



Swan  
representation


Swan represents transcript models as a series of splice sites and intronic or exonic connections between them

- Assign each unique splice site to a node
- Nodes are colored by their role in transcript model (TSS, TES)




Swan represents transcript models as a series of splice sites and intronic or exonic connections between them

- Assign each unique splice site to a node
- Nodes are colored by their role in transcript model (TSS, TES)
- Connect nodes (splice sites) by intronic or exonic edges




Swan represents transcript models as a series of splice sites and intronic or exonic connections between them

- Assign each unique splice site to a node
- Nodes are colored by their role in transcript model (TSS, TES)
- Connect nodes (splice sites) by intronic or exonic edges



Swan represents transcript models as a series of splice sites and intronic or exonic connections between them


- Assign each unique splice site to a node
- Nodes are colored by their role in transcript model (TSS, TES)
- Connect nodes (splice sites) by intronic or exonic edges





# Swan represents transcript models as a series of splice sites and intronic or exonic connections between them

- Assign each unique splice site to a node
- Nodes are colored by their role in transcript model (TSS, TES)
- Connect nodes (splice sites) by intronic or exonic edges



# Swan visualizations: Gene summary plots and transcript path plots


**nodes:** TSS TES

internal


**edges:** exon intron

**dashed edges:** novel  
intron or exon

**outlined nodes:** novel  
splice site, TSS, or TES




*Tfdpi1* Gene Summary Plot




*Tfdpi1* Transcript Path Plot

# Data and data analysis overview

- Publicly available PacBio RNA-seq data
- C2C12 skeletal muscle cell line model of myogenesis
  - 0hr - undifferentiated, mononucleated myoblasts
  - 72 hr - differentiated, multinucleated myotubes
- Found isoform-switching, TSS-switching, and TES-switching genes with Swan




C2C12 0hr myoblasts




C2C12 72hr myotubes


# Isoform-switching gene *Chp1* report




# Isoform-switching gene *Chp1* report




# TSS-switching gene *Tnnt2* report



# TSS-switching gene *Tnnt2* report




# TES-switching gene *Tpm1* report







# TES-switching gene *Tpm1* report




# Swan's representation of TSSs, TESs, and isoforms are compatible with Scanpy plotting tools



Isoform expression




TSS expression




TES expression

# Swan can be used to analyze single-cell transcriptomes





Clusters and cell type labels of 464 differentiating C2C12 skeletal muscle cells and nuclei

# Swan can be used to analyze single-cell transcriptomes



*Tpm2* gene report in differentiating C2C12 skeletal muscle


# Swan can be used to analyze single-cell transcriptomes




Isoform and TSS expression plots in differentiating C2C12 skeletal muscle

# Your turn!

- Publicly-available bulk PacBio RNA-seq data from 0hr (myoblast) and 72hr differentiated (myotube) C2C12 skeletal muscle
- Reads cleaned with TranscriptClean and annotated with TALON
- DEGs, DETs, and isoform-switching, TSS-switching, and TES-switching genes called with Swan
- Novel intron retention and exon skipping events called with Swan



C2C12 0hr myoblasts



C2C12 72hr myotubes

# Acknowledgements

## Mortazavi Lab:

Elisabeth Rebboah

Gaby Balderrama-Gutierrez

Dana Wyman, PhD (former)

Katherine Williams

Heidi Liang

Cassandra McGill (former)

Klebea Carvalho

Christina Wilcox


Narges Rezaie

Jasmine Sakr

Angelique Cortez


Isaryhia Rodriguez

Hasan Celik



## Wold Lab (Caltech):

Diane Trout



Center for Complex Biological Systems  
University of California, Irvine

NHGRI UM1 HG009443