


Single Cell Sample Preparation




Diana Burkart-Waco, PhD

dburkart@ucdavis.edu

SRA at UC Davis DNATECH Core

<https://dnatech.genomecenter.ucdavis.edu/>

General Workflow



Cell isolation

- Cell isolation guides available at:

<https://www.support.10xgenomics.com/single-cell-gene-expression/sample-prep/>



Sample Prep

- Demonstrated Protocol (14 documents)
 - Single Cell Gene Expression Demonstrated Protocol Compatibility Table
 - Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols
 - Methanol Fixation of Cells for Single Cell RNA Sequencing
 - Isolation of Nuclei for Single Cell RNA Sequencing
 - Single Cell Protocols - Cell Preparation Guide
 - Enrichment of CD3+ T Cells from Dissociated Tissues for Single Cell RNA Sequencing and Immune Repertoire Profiling
 - Tumor Dissociation for Single Cell RNA Sequencing
 - Thawing Dissociated Tumor Cells for Single Cell RNA Sequencing
 - Single Cell Suspensions from Cultured Cell Lines for Single Cell RNA Sequencing
 - Removal of Dead Cells from Single Cell Suspensions for Single Cell RNA Sequencing
 - Moss Protoplast Suspension for Single Cell RNA Sequencing
 - Fresh Frozen Human-Mouse Cell Line Mixtures for Single Cell RNA Sequencing
 - Fresh Frozen Human Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing
 - Dissociation of Mouse Embryonic Neural Tissue for Single Cell RNA Sequencing

- Cell isolation → biggest source technical variation.
- Dissociation and preparation depends on cell type.

Tissue Dissociation

- Many ways to dissociate tissue:
 - Mechanical
 - Enzymatic
 - Automated blending
 - Microfluidics devices
- Considerations:
 - Speed
 - Consistency in results
 - Good representation of all cell types




Arabidopsis
protoplasts

Same method,
different days.

Factors Influencing Dissociation

DNA
TECH



Unfortunately, SC sample optimization is best achieved through trial and error...


Single cell vs. nuclei

- Single cell captures more transcripts, but is a harder protocol.
- When to use nuclei:
 - Cells cannot be harvested intact or viable (e.g. adipocytes, neurons).
 - Cells are too big for capture (e.g. cardiomyocytes).
 - Tissue frozen.

Tissue dissection



Cell lysis



Nuclei purification



Detergent
Grinding

Flow cytometry
Density Gradient

Many protocols available...



- Start here:

Customer Developed Protocol

'Frankenstein' protocol for nuclei isolation from
fresh and frozen tissue



Contributed by:

Luciano Martelotto, Ph.D., Melbourne, Centre for Cancer
Comprehensive Cancer Centre



www.collaslab.com

Isolation of Nuclei from Somatic Cells

1. HeLa Cells, 293T Cells, NT2 Cells



Cell preparation

- harvest cells from flasks as per standard protocol
- spin cells in 50 ml conical tube at 1500 rpm for 10 min at RT
- resuspend cells in 30 ml PBS; take a 50 μ l sample to determine concentration



10x Genomics®

Sample Preparation Demonstrated Protocols


Isolation of Nuclei for Single Cell RNA
Sequencing

Factors influencing success


- Cell debris / dead cells
 - Clog microfluidics, free floating RNA → noise.
- Aggregates
 - No longer single cell data.
 - Clog microfluidics.
- Buffer
 - Inhibit polymerase → decrease library complexity.
- Storage conditions

But the most important factor is **cell / nuclei counting!**


Workflow



Tissue aggregates



Tissue clump




Single nucleus

DAPI stain
60X
No scale bar

Recommended treatment: optimize tissue dissociation

Cell debris

- Mouse DRG.



Clean DRG sample



Noisy sample

Recommended treatment: filtration, centrifuge, add blocking agent.

Cell debris II

- Lettuce nuclei prep.






Photo credit: Mohan Prem Anand Marimuthu

Recommended treatment: modify dissociation times / detergent concentrations, change density gradient.

Dead cells (+debris?)




Recommended treatment: dead cell removal.

<https://www.miltenyibiotec.com/US-en/products/dead-cell-removal-kit.html>

Doublets

- Non-single cell clumps.
- Integrated into droplets and cannot* distinguish from single cells.




Bad for any
single cell
experiment.




Recommended treatment: filtration, change [blocking agent]

Workflow



Automated counting



- Countess II.



- Pros (+):
 - Fast.
 - Live/dead cell counts.
 - Cell size estimates.
- Cons (-):
 - Cell size limits (4-30um).
 - Doesn't do well with odd cell shapes, debris, nuclei.

Manual counting

- Hemocytometer.





- Pros (+):
 - Reliable cell counts.
 - Count small cells.
 - Visualize cells vs. debris.
- Cons (-):
 - Slow.

Nuclei QC

- Use fluorescence.
 - Countess II FL, Devonix CellDrop, Luna FL.
- Stains
 - Trypan blue → 100% dead (not great).
 - Ethidium homodimer stain (good for excluding debris).
 - DAPI, PI other options




60X: DAPI



Credit: 10X

Cell sorting





- FACS





- Pros (+):
 - Sort live / dead.
 - Enrich cells of interest.
 - Determine whether correct cells isolated (qPCR works too).
- Cons (-):
 - Overestimates cell counts (~2x).

<https://flowcytometry-embl.de/cell-sorting/>

The good...





The bad...





Challenges Cell Recovery

- 10X is optimized for model cell types – human PBMCs, etc.
 - Round, easy to count, size well below size of microfluidics.
 - These yield consistent results. 10K target. 8K-10K recovered.



- But most experiments outside of culture don't look like this.



Variable size,
shape, and
viability

How maximize output?


- Carefully craft **experimental design** and **sample prep** .
- QC **cells** before real sample set-up.
- Concentration: aim for the median.
 - **700-1,200 cells per µl optimal.**
 - Too high → dilute.
 - Too low → tough one (concentration impacts yield, pooling replicates suboptimal).
 - Count in replicates (at least n=2).
- Viability: **70% minimum.**
 - Nuclei and methanol fixed cells (0%).
- Treat cells gently.
 - Wide bore pipette tips, keep cells at preferred temp, work quickly.

10X Single Cell Studies - DNATECH



Assay summary (organism)

- 78 projects since July 2017 (primarily 3' GEX).



Sample
distribution by
organism

Aurelia aurita
Bovine
Chicken
Fruit fly
Human
Human+Mouse
Hydra
Mouse
Rhesus macaque
Zebrafish

SCIENTIFIC REPORTS



OPEN

Molecular profiling of resident and infiltrating mononuclear phagocytes during rapid adult retinal degeneration using single-cell RNA sequencing

Received: 15 October 2018

Accepted: 27 February 2019

Published online: 19 March 2019

Kaitryn E. Ronning¹, Sarah J. Karlen², Eric B. Miller¹ & Marie E. Burns^{1,2,3}

Vitamin E deficiency - DRG Neurons



CellPress

Sneak Peek

A PREVIEW OF PAPERS UNDER REVIEW

Single-Cell RNA-Seq Reveals Profound Alterations in Mechanosensitive Dorsal Root Ganglion Neurons with Vitamin E Deficiency

iScience

50 Pages • Posted: 5 Aug 2019 • Sneak Peek Status: **Review Complete**



Carrie Finno

University of California, Davis, School of Veterinary Medicine, Department of Population Health and Reproduction

Janel Peterson

University of California, Davis, School of Veterinary Medicine, Department of Population Health and Reproduction

'564' – challenges and solutions



- DRG fragile, difficult to count, odd sizes (clog 10X chip?).
- Count manually → better counts.
- Accepted some degree of induced technical variation
 - Prepped samples over three days.
 - Improved sample quality and reproducibility, introduced prep day variation.
 - As Matt mentioned mitigate with experimental design.

Important resources

- 10X Genomics
 - <https://support.10xgenomics.com/single-cell-gene-expression>
 - Paul Scott, Sales Executive (paul.scott@10xgenomics.com)
- UC Davis Flow Cytometry
 - http://www.ucdmc.ucdavis.edu/pathology/research/research_labs/flow_cytometry/index.html
 - Bridget McLaughlin (Technical Director)
- UC Davis DNA Technology Core
 - <http://dnatech.genomecenter.ucdavis.edu/single-cell-analyses/>
- UC Davis Bioinformatics Core
 - <https://bioinformatics.ucdavis.edu/>


Single Cell Expression Profiling & Genomics (10X Genomics, Fluidigm & Plate-based scRNA-Seq)

10X Genomics Single Cell Sequencing

Please schedule any 10X single cell experiment at least a week in advance. We highly advise a consultation prior to experiment scheduling.

Update March 2019: MULTI-seq reagents promise to enable the labeling and pooling of up to 96 cell suspension samples. Reagents are in stock. Please inquire!

The 10X Genomics Single Cell system is the **single-cell expression profiling** platform enabling the analysis of large cell numbers at the highest capture efficiency (of up to 65%). The technology allows for high-throughput single cell transcriptomics of many different cell types as well as **single-nuclei expression profiling**. The flexible workflow encapsulates 500 to 20,000 cells or nuclei per library together with micro-beads into nano-droplets. Each bead is loaded with adapters containing one of 750,000 different barcodes for the single cell RNA-seq library prep. In contrast to other protocols (e.g. Drop-Seq) the 10X controller is capable of loading "all" droplets with micro-beads, enabling single-Poisson distribution loading and thus high capture efficiencies (in contrast to double-Poisson loading of other protocols). The single-cell encapsulating process is significantly faster compared to *inDrop* or *Drop-Seq*. Up to eight samples can be processed per batch within minutes. The resulting data can be analyzed with the free *Cell Ranger* and *Loupe Cell Browser* software. In addition the **Bioinformatics Core** has developed a custom single-cell data analysis pipeline for 10X data.



Cost 10X Reagents

Product Code	Product	List Price	Discount
1000120	Chromium Next GEM Chip G Single Cell Kit, 48 rxns	USD 1,480.00	
1000121	Chromium Next GEM Single Cell 3' GEM, Library & Gel Bead Kit v3.1, 16 rxns	USD 21,600.00	

Quote Created Date

3/10/2020

- Chips are single use (n=1 to n=8 samples) → \$247 EACH
- Reagents (n=1) → \$1,350 EACH.
- Total (n=1, reagents only): \$1,597.
 - This doesn't include you or your cores labor!
- Sequencing? Bioinformatics?
 - Depends on many factors.
 - We charge \$1,000 for 350M reads (~10K cells with 35K reads per cell).

Acknowledgements

- UC Davis DNA Technology Core team.
- UC Davis Bioinformatics Core team.
- 10X Genomics.
 - Adam Bemis (Field Applications Scientist)
 - Paul Scott (Sales Executive)
 - Nicole Abreu (Technology Advisor)

