



PACIFIC  
BIOSCIENCES®




# Assembling the drosophila genome with IPA and HiFi data

Zev Kronenberg

---

# AGENDA

- Introduction
- Setup PacBio software
- Learn about command line interface
- Run an IPA assembly




# Today's dataset

Drosophila

# DROSOPHILA DATASET INFO

## 19 Kb dataset

- Processing of PacBio data
  - Circular consensus algorithm was done with SMRT Link
  - Data was subsampled down to 38x depth of coverage (DOC)
  
- Short read data
  - Standard procedures
  - Both parental strains were sequenced (70-90x DOC)
  - Utility
    - Trio binning
    - Phasing evaluation



# SETUP CONDA – A SOFTWARE MANAGEMENT SYSTEM

- Google: `conda install`
  - <https://docs.conda.io/projects/conda/en/latest/user-guide/install/linux.html>
- Follow link for download (linux x86)
- Follow install directions, **setup BASH shell!**
- **You will likely need to source your .bashrc if conda isn't in your path after setup**

- Why does PacBio use Conda?
  - Central repository of many bioinformatic tools
  - We can distribute binary code
  - We can post updates quickly
  - Many more...

## SETUP IPA



PACBIO® & BIOCONDA®

- Google: `pacbio bioconda`
  - <https://github.com/PacificBiosciences/pbbioconda>
  - Familiarize yourself with available command line tools available for download.
  - Go to IPA bioconda wiki:
  - <https://github.com/PacificBiosciences/pbbioconda/wiki/Improved-Phased-Assembler>



## SETUP IPA



PACBIO®



BIOCONDA®

```
conda create -n ipa -c  
bioconda -c conda-forge -c  
defaults conda activate ipa  
conda install pbipa
```



## UC DAVIS SETUP



PACBIO®

& BIOCONDA®

```
eval "$(/share/biocore/shunter/2020-07-15-IPA-tests/conda/bin/conda  
shell.bash hook)"
```





# COMMAND LINE



PACBIO®



BIOCONDA®

```
(ipa) zevk@tadpole:~$ ipa -h
usage: ipa [-h] [--version] {local,dist,validate} ...

Improved Phased Assembly tool for HiFi reads.

optional arguments:
  -h, --help            show this help message and exit
  --version              show program's version number and exit

subcommands:
  One of these must follow the options listed above and may be followed by sub-command specific options.

  {local,dist,validate}
                        sub-command help
  local                Run IPA on your local machine.
  dist                 Distribute IPA jobs to your cluster.
  validate              Check dependencies.

Try "ipa local --help".
Or "ipa validate" to validate dependencies.
https://github.com/PacificBiosciences/pbbioconda/wiki/Improved-Phased-Assembler
```

# KNOW YOUR VERSIONS



PACBIO®



BIOCONDA®

- We commonly update IPA
- Before running assembly update IPA
- Keep track of your versions
- Avoid updating mid-assembly

```
(ipa) zevk@tadpole:~$ ipa validate
INFO: /home/zevk/anaconda3/envs/ipa/bin/ipa validate
Checking dependencies ...
/home/zevk/anaconda3/envs/ipa/bin/python3
/home/zevk/anaconda3/envs/ipa/bin/ipa2-task
/home/zevk/anaconda3/envs/ipa/bin/falconc
/home/zevk/anaconda3/envs/ipa/bin/nighthawk
/home/zevk/anaconda3/envs/ipa/bin/pancake
/home/zevk/anaconda3/envs/ipa/bin/pblayout
/home/zevk/anaconda3/envs/ipa/bin/racon
/home/zevk/anaconda3/envs/ipa/bin/samtools
snakemake version=5.20.1
Machine name: 'Linux'
ipa2-task 0.2.0 (commit 33ccb062c1db781cd9aa10e4341c670430b1e575)
falconc version=1.5.1+git.895d7f33113c17b399428ff45dce127f7aa635ef, nim-version=1.2.0
Nighthawk 0.1.0 (commit df65ce5*)
pancake 0.1.0 (commit 3a4146f*)
pblayout 0.1.0 (commit 5257a1a*)
racon version=v1.4.13
samtools 1.9
Using htlib 1.9
```

# RUN IPA ON THE CLUSTER



PACBIO®



BIOCONDA®

```
ipa dist -i ../hifi_long_read_data/ELF_19kb.m64001_190914_015449.Q20.38X.fasta \  
--nthreads 24 --njobs 30 --cluster-args 'sbatch -J zev-ipa.{rule} -t 45 \  
-c {params.num_threads} -e stderr -o stdout --get-user-env \  
--chdir pacbio_2020_data_drosophila/hifi_long_read_diploid_ipa_assembly_cluster'
```

## RUN IPA LOCALLY



```
ipa local --nthreads 48 --njobs 2 -i ELF_19kb.m64001_190914_015449.Q20.38X.fasta
```



PACBIO®



BIOCONDA®

# STAGES OF IPA

```
drwxrwsr-x 2 zevk genome_workshop 7 Jul 14 14:17 01-generate_config
drwxrwsr-x 2 zevk genome_workshop 10 Jul 14 14:17 02-build_db
drwxrwsr-x 3 zevk genome_workshop 3 Jul 14 14:18 03-ovl_asym_prepare
drwxrwsr-x 5 zevk genome_workshop 5 Jul 14 14:21 04-ovl_asym_run
drwxrwsr-x 2 zevk genome_workshop 8 Jul 14 14:24 05-ovl_asym_merge
drwxrwsr-x 3 zevk genome_workshop 3 Jul 14 14:24 06-phasing_prepare
drwxrwsr-x 8 zevk genome_workshop 8 Jul 14 14:43 07-phasing_run
drwxrwsr-x 2 zevk genome_workshop 10 Jul 14 14:49 08-phasing_merge
drwxrwsr-x 2 zevk genome_workshop 111 Jul 14 14:52 09-ovl_filter
drwxrwsr-x 2 zevk genome_workshop 36 Jul 14 14:57 10-assemble
drwxrwsr-x 3 zevk genome_workshop 3 Jul 14 14:57 11-polish_prepare
drwxrwsr-x 5 zevk genome_workshop 5 Jul 14 15:04 12-polish_run
drwxrwsr-x 2 zevk genome_workshop 7 Jul 14 15:10 13-polish_merge
drwxrwsr-x 2 zevk genome_workshop 4 Jul 14 15:10 14-final
-rw-rw-r-- 1 zevk genome_workshop 210 Jul 14 14:17 config.json
-rw-rw-r-- 1 zevk genome_workshop 140 Jul 14 14:17 config.yaml
-rw-rw-r-- 1 zevk genome_workshop 124 Jul 14 14:17 input.fofn
-rw-rw-r-- 1 zevk genome_workshop 307 Jul 14 14:17 ipa.log
drwxrwsr-x 2 zevk genome_workshop 2 Jul 14 14:17 qsub_log
```

# GETTING ASM STATS



PACBIO®



BIOCONDA®

```
module load assembly_stats/1.0.1
```

```
assembly-stats -t final.p_ctg.fasta  
final.a_ctg.fasta | column -t
```

filename	total_length	number	mean_length	longest	shortest	N_count	Gaps	N50	N50n	N70	N70n	N90	N90n
../hifi_long_read_diploid_ipa_assembly/RUN/14-final/final.p_ctg.fasta	232197789	208	1116335.52	23464522	50222	0	0	7970785	9	2631370	18	461442	54
../hifi_long_read_diploid_ipa_assembly/RUN/14-final/final.a_ctg.fasta	36577756	287	127448.63	9995912	4424	0	0	1540414	5	649690	11	28408	128

## LOOKING AT YOUR ASM



PACBIO®



BIOCONDA®

Make GFA

```
pblayout subgraph --tp all_ctg.tp sg_edges_list tp_graph
```

Remove

Double


Edges

```
perl -lane '$_ =~ s/:B|:E//g; print' tp_graph.csv >  
tp_graph_single.csv ; perl -lane '$_ =~ s/:B|:E//g;  
print' tp_graph.gfa > tp_graph_single.gfa
```

# LOOKING AT YOUR ASM




PACBIO® & BIOCONDA®



## Bandage

a Bioinformatics Application for Navigating *De novo* Assembly Graphs Easily

Bandage is a program for visualising *de novo* assembly graphs. By displaying connections which are not present in the contigs file, Bandage opens up new possibilities for analysing *de novo* assemblies.



Download Mac

Download Linux

Download Windows

View project on GitHub

For installation







PACBIO®

[www.pacb.com](http://www.pacb.com)

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.