



PACIFIC
BIOSCIENCES®

Improved Phased Assembly using HiFi Data

Ivan Savić, Ph.D., PacBio Assembly Tech Lead

Zev Kronenberg, Christopher Dunn, Derek Barnett, Sarah Kingan, James Drake, Jonas Korlach

UC Davis Workshop, 2020



@IvanSovic

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. All other trademarks are the sole property of their respective owners.

IMPROVED PHASED ASSEMBLY USING HIFI DATA

James Drake

Derek Barnett

Ivan Sović

Zev Kronenberg


Christopher Dunn



Sarah Kingan





Jonas Korlach



WHAT ARE HIFI READS?

- **They are long**
 - Up to 25 kb
- **They are accurate**
 - Long reads with $\geq Q20$ (99%) accuracy
- **They have single-molecule resolution**
 - Sequence DNA or RNA
- **They have little bias**
 - No DNA amplification, least GC content and sequence complexity bias



HiFi READ
(>99% Accuracy)

[illegible]

19,820 bp HiFi read, predicted QV: 33
19,812 bp correct, 8 errors
 99.96% accurate (QV34)


98.16% accuracy

98.16% accuracy

HIFI READS FOR IMPROVED ASSEMBLY

Contiguity

- Resolve Repetitive Regions
- High Contig N50



Correctness


- Base QV
- Phasing accuracy

AGTTTCGATAGA

AGTT-CGAA**AGA**

Completeness

- Gene Space
- Repetitive Regions




Compute

- CPU / Wall Time
- RAM
- Disk Storage



ASSEMBLY METHODS FOR HIFI READS


FALCON-Unzip



hifiasm


 **Heng Li**
@lh3lh3

@ChengChhy from my group developed hifiiasm, a new #Hifireads assembler that preserves local phasing. Human asm in 1/2 day. Tested on non-humans. Comparable asm quality to others. Features in plan: consistent primary asm & global phasing. Feedback welcomed!

chhylp123/hifiiasm
Hifiiasm: a haplotype-resolved assembler for accurate Hifi
reads - chhylp123/hifiiasm
 [github.com](https://github.com/chhylp123/hifiiasm)

8:58 AM · Jan 14, 2020 · Twitter Web App

HiCanu



Flye

Peregrine




<https://github.com/cschin/Peregrine>
<https://github.com/chhyip123/hifiasm>
<https://github.com/fenderglass/Flye>
<https://www.biorxiv.org/content/10.1101/2020.03.14.992248v3>
<https://www.pacb.com/proceedings/high-quality-human-genomes-achieved-through-hifi-sequence-data-and-falcon-unzip-assembly/>

IMPROVED AND PHASED ASSEMBLY (IPA)




IMPROVED AND PHASED ASSEMBLY (IPA)



Goals:

1. Fast assembly and quick turnaround time
2. High contiguity
3. Fully phased haplotigs
4. High per-base quality of polished assemblies
5. Ease of use

IMPROVED AND PHASED ASSEMBLY (IPA)



—Work in progress


- Currently in Beta
- Rapidly being updated

IPA METHODS

IPA WORKFLOW




IPA WORKFLOW



—Sequence Database


IPA WORKFLOW



— Sequence Database

- Converting one or more input files into a unified database format
- **SeqDB** – database of all input reads, for fast random access
- **SeedDB** – database of seeds (e.g. minimizers) precomputed from the DB


IPA WORKFLOW



— Sequence Database


- Supported formats: FASTA, FASTQ, BAM, XML and FOFN (including gzipped FASTA and FASTQ)
- Compression
- Arbitrary method for seed generation in SeedDB
 - Minimizers
 - Full set of dense k-mers
 - **Spaced seeds**
 - Other approaches are trivial to add
- Other features:
 - Fetching sequences in original or homopolymer compressed space
 - Converting headers to IDs

IPA WORKFLOW



—Pancake


IPA WORKFLOW



—Pancake


- New overlapper
- Extremely fast and accurate
- Can overlap a 30x NA19240 (18kb) dataset in **20 CPU hrs**

IPA WORKFLOW




—Pancake

- Runs on a pair of blocks from the SeqDB (query and target)
- Algorithm:
 - For each read in the query block collect all seed hits in the target block
 - Sort and bin the seed hits in narrow diagonal bins
 - Initial seeding of potential local alignments
 - For each candidate diagonal, perform fast alignment computation and alignment extension to form dovetail overlaps
 - Filter low quality overlaps




IPA WORKFLOW




—Pancake

- Spaced seeds (minimizers)
 - Novel adaptation in combination with minimizers
 - Seeds are constructed by skipping zero or more bases after every inclusive base
 - Efficient to compute
 - Seeds cover larger regions
 - Increases specificity
 - Minimizer approach applied on spaced seeds




- By default, spacing of 1 is used

IPA WORKFLOW



— Nighthawk


IPA WORKFLOW



— Nighthawk

- New phasing tool
- Novel approach based on the de Bruijn graph!
- Works directly on overlap piles!

IPA WORKFLOW




— Nighthawk


— Idea:

- Discover and remove overlaps between reads coming from different haplotypes
- **Phasing before layout – unlike FALCON-Unzip**
- Goal: Natural phase separation at the layout stage

Assembly string graph - HiFi read overlaps



IPA WORKFLOW




— Nighthawk

— Algorithm:

- Builds a de Bruijn graph for each overlap pile
- Analyzes the bubbles
- Computes a Read Similarity Score for each pair of reads
- Phases bubbles in the de Bruijn graph
- Performs transitive inference
- Finally – filters cross-phase overlaps


A.

Overlap reads (overlap pile - m4)




B.

Build read-colored de Bruijn graph




C.

Score all read overlaps - per pile




D.

Filter overlaps and run std. assembly




IPA WORKFLOW




— Nighthawk

- Visualization of a de Bruijn graph for a pile of reads of a *Drosophila* HiFi dataset
 - $k = 23$
 - Red dots – heads of heterozygous bubbles



IPA WORKFLOW



— Nighthawk

- More info on Nighthawk in our blog post:
 - <https://www.pacb.com/blog/direct-phased-genome-assembly-using-nighthawk-on-hifi-reads/>

Direct Phased Genome Assembly Using Nighthawk on HiFi Reads

Monday, January 13, 2020

By Zev Kronenberg, Senior Engineer of Bioinformatics at PacBio

Since the introduction of HiFi reads the community has embraced these long and highly accurate reads for human genome assembly and paralog resolution [1-5]. At PacBio, the assembly team (Figure 1) is working to build on the accuracy of HiFi data for direct phasing during assembly.






Figure 1. The PacBio assembly team. From left to right, James Drake, Zev Kronenberg (@ZevKronenberg), Derek Barnett (@DerekWBarnett), Chris Dunn, and Ivan Sović (@IvanSovic)

IPA WORKFLOW



— Chimera and repeat filtering


IPA WORKFLOW




— Chimera and repeat filtering

- Small fraction of HiFi reads are molecular chimeras
- Filtering improves contiguity and reduces misassemblies


Visualization of chimeric joins in an assembly graph



Overlap pile demonstrating chimera detection



IPA WORKFLOW



—Layout


IPA WORKFLOW




—Layout

- String graph based
- Polyploid aware
- Primary and associate contig sets
- Phase-aware read tracking
 - Reads assigned to contigs based on phased overlaps
 - Important for polishing

String graph for a Drosophila HiFi dataset




IPA WORKFLOW



—Polishing

IPA WORKFLOW




— Polishing

- Consumes read-to-contig assignment
- Phase-aware
- Assignment-based mapping using **Pbmm2**

— Racon

- Possibility of GPU acceleration



<https://github.com/isovic/racon>

IPA WORKFLOW

— Full workflow with phasing and polishing



IPA WORKFLOW

— Haploid workflow – phasing can optionally be switched on/off



IPA WORKFLOW

- Polishing can optionally be switched on/off
- Fast draft assembly

Phased workflow



Haploid workflow




RESULTS

RESULTS: HUMAN ASSEMBLY IS VERY FAST

HPRC HG002 34x Dataset – Haploid workflow without polishing


	FALCON	IPA
N50 [Mbp]	31.37	38.81
Max length [Mbp]	110.15	110.72
Total length [Gbp]	2.96	3.06
CPU time [h]	2767	46
		60x Faster!



RESULTS: LONG PHASE BLOCKS IN HUMAN, HIGH BASE QV

HPRC HG002 34x Dataset – Phased workflow with polishing

	FALCON-Unzip		IPA (Phased)	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	31.40	0.191	33.75	0.352
Max length [Mbp]	110.12	1.62	110.94	2.30
Total length [Gbp]	2.95	1.99	3.02	1.85
CPU time [h]	5102		590	
			8.64x Faster!	



RESULTS: LONG PHASE BLOCKS IN HUMAN, HIGH BASE QV

HPRC HG002 34x Dataset – Phased workflow with polishing

	FALCON-Unzip		IPA (Phased)	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	31.40	0.191	33.75	0.352
Max length [Mbp]	110.12	1.62	110.94	2.30
Total length [Gbp]	2.95	1.99	3.02	1.85
CPU time [h]	5102		590	
			8.64x Faster!	


- Primary contig pile is slightly larger than expected haploid genome size
- Fully phased regions of the graph can appear as separate graph components

RESULTS: LONG PHASE BLOCKS IN HUMAN, HIGH BASE QV


Purging “duplicate” haplotigs from the primary contig set

- Common expectation when phasing contigs


(a) FALCON




(b)



(c) FALCON-Unzip




Example of a phased assembly graph




RESULTS: LONG PHASE BLOCKS IN HUMAN, HIGH BASE QV

Purging “duplicate” haplotigs from the primary contig set

Fully phased graph component




Artifacts in the graph - spurs



RESULTS: LONG PHASE BLOCKS IN HUMAN, HIGH BASE QV

Purging “duplicate” haplotigs from the primary contig set

- Happens to all current assembly tools
- Remedy – publicly available tool “purge_dups”



The screenshot shows the GitHub repository page for `dfguan/purge_dups`. The repository has 32 issues, 1 pull request, 1 project, 1 wiki, 1 security, and 1 insight. The main branch is `master` with 6 branches and 3 tags. The repository contains the following files:

File	Commit	Time
<code>scripts</code>	<code>v124-master: a bug in run_purge_dups.py, input filenames with sam...</code>	23 days ago
<code>src</code>	<code>v123-master: print wrong peak numbers (not affecting results)</code>	3 months ago
<code>LICENSE</code>	<code>Update LICENSE</code>	12 months ago
<code>README.md</code>	<code>v124-master: update readme</code>	12 days ago
<code>purge_duppipeline.png</code>	<code>v3: update purge_duppipeline</code>	14 months ago

The `README.md` file is titled **Purge_Dups** and describes the tool as a haplotypic duplication identification tool. The repository also has 77 commits and 1 release (v1.0.1, Latest, on 28 Feb).

About
haplotypic duplication identification tool
Readme
MIT License

Releases 3
v1.0.1 (Latest) on 28 Feb
+ 2 releases

Languages
C 86.3% Python 11.3% Shell 2.0% Makefile 0.4%

RESULTS: GREAT HAPLOTIG SEPARATION WITH PURGE DUPS

HPRC HG002 34x Dataset – Phased workflow with polishing

	FALCON-Unzip + Purge dups		IPA (Phased) + Purge dups	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	33.25	0.195	34.48	0.353
Max length [Mbp]	110.12	1.62	110.94	4.12
Total length [Gbp]	2.87	1.98	2.88	1.94
Base QV	50.6	49.9	50.6	50.2
Phase accuracy	0.706	0.997	0.720	0.980
BUSCO of primary	C:95.1% S:94.2%,D:0.9%		C:95.2% S:94.2%,D:1.0%	


RESULTS: GREAT HAPLOTIG SEPARATION WITH PURGE DUPS

HPRC HG002 34x Dataset – Phased workflow with polishing

	FALCON-Unzip + Purge dups		IPA (Phased) + Purge dups	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	33.25	0.195	34.48	0.353
Max length [Mbp]	110.12	1.62	110.94	4.12
Total length [Gbp]	2.87	1.98	2.88	1.94
Base QV	50.6	49.9	50.6	50.2
Phase accuracy	0.706	0.997	0.720	0.980
BUSCO of primary	C:95.1% S:94.2%,D:0.9%		C:95.2% S:94.2%,D:1.0%	

RESULTS: HIGHLY ACCURATE CONTIG ASSEMBLY

Atlantic Bluefin Tuna – Phased workflow




	IPA (Phased)		IPA + purge_dups	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	9.34	3.70	13.80	3.83
Max length [Mbp]	39.38	13.95	39.38	19.49
Total length [Gbp]	1.26	0.280	0.791	0.744
BUSCO of primary	C:97.4% S:44.6%,D:52.8%		C:97.7% S:95.1%,D:2.6%	

* Credit for the ABFT data, reference assembly and scaffolding: Paul Peluso¹, Greg Concepcion¹, Jay Ghurye², Nathan Truelove³ and Barbara Block⁴

¹ Pacific Biosciences, ² Dovetail Genomics, ³ Monterey Bay Aquarium Research Institute,


⁴ Hopkins Marine Station




* Cabanettes F, Klopp C. (2018) D-GENIES: dot plot large genomes in an interactive, efficient and simple way. PeerJ 6:e4958 <https://doi.org/10.7717/peerj.4958>

RESULTS: HIGH PHASE ACCURACY

Drosophila melanogaster F1 – Phased and polished
















	Hifiasm + purge_dups		IPA + purge_dups	
	primary	haplotigs	primary	haplotigs
N50 [Mbp]	22.55	1.28	13.49	2.42
Max length [Mbp]	28.13	6.81	23.47	12.48
Total length [Mbp]	160.19	149.87	134.19	115.26
Base QV	48.1	47.4	47.97	46.87
Phase accuracy	0.788	0.998	0.826	0.999
BUSCO of primary	C:98.5% S:98.0%,D:0.5%		C:98.7% S:98.2%,D:0.5%	



RESULTS: BUG GENOMES

Testing on real-world samples - butterflies, moths & mosquitoes
Darwin Tree of Life, Sanger



RESULTS: BUG GENOMES





ALL SAMPLES: IPA vs. FALCON

Species	Primary										Haplotigs									
	# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV		# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA
ilAutoGam1	85	79	10.97	12.01	381	368	99.0%	99.2%	47.1	47.1	378	1613	6.98	9.09	330	379	93.8%	97.9%	46.3	44.9
ilCosmTra1	2173	1924	0.79	0.88	872	862	97.9%	97.5%	44.4	45.3	3932	3505	0.27	0.32	717	742	85.6%	89.8%	44.1	43.3
ilCranLig1	220	267	7.32	3.50	438	436	99.0%	98.9%	46.1	46.8	4218	1241	0.07	0.42	246	266	57.3%	67.7%	39.1	45.6
ilEndoFla1	623	489	1.46	1.92	492	489	98.6%	99.2%	46.5	45.9	2110	1587	0.29	0.56	375	418	82.7%	89.0%	45.9	46.3
ilLymaMon1	251	301	10.31	5.77	917	912	99.2%	99.2%	46.5	46.8	5525	3162	0.17	0.47	610	633	64.6%	66.4%	41.6	43.8
ilNoctFim1	382	307	3.12	3.75	576	572	98.2%	98.9%	47.6	46.3	1493	1152	0.86	1.62	514	545	83.9%	94.9%	45.1	48.1
ilNotoUdd1	1128	809	1.52	2.07	829	814	98.9%	98.3%	45.5	46.8	4115	2791	0.26	0.69	629	780	75.3%	94.0%	42.8	45.6
ilParaStr1	218	129	4.42	6.92	480	481	99.4%	99.4%	48.0	47.4	1166	1309	0.64	1.78	413	431	84.3%	84.7%	47.1	46.9
ilRecuLeu1	1268	1247	1.09	0.98	748	746	98.5%	98.2%	44.4	45.1	4890	4029	0.12	0.23	416	534	54.1%	70.4%	42.0	43.4
ilThyaBat1	186	88	3.29	7.15	316	316	98.5%	98.9%	46.2	47.0	650	1073	0.89	2.77	296	327	91.6%	96.0%	46.9	45.3
ilVaneAta1	80	48	10.13	12.12	368	368	99.1%	99.3%	47.9	48.3	2185	1444	0.14	4.80	205	369	60.1%	97.7%	45.3	46.9
idAnopAqu88	260	79	18.22	15.05	188	181	99.2%	97.8%	49.5	50.6	1056	1833	0.11	4.02	90	211	43.2%	94.2%	48.9	43.8
idAnopCol22p13	272	189	5.48	5.01	260	260	99.4%	99.2%	49.4	49.5	1112	2474	0.12	0.23	77	148	27.1%	52.5%	43.9	40.3

RESULTS: BUG GENOMES

ALL SAMPLES: IPA vs. FALCON

Species	Primary										Haplotigs									
	# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV		# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA
ilAutoGam1	85	79	10.97	12.01	381	368	99.0%	99.2%	47.1	47.1	378	1613	6.98	9.09	330	379	93.8%	97.9%	46.3	44.9
ilCosmTra1	2173	1924	0.79	0.88	872	862	97.9%	97.5%	44.4	45.3	3932	3505	0.27	0.32	717	742	85.6%	89.8%	44.1	43.3
ilCranLig1	220	267	7.32	3.50	438	436	99.0%	98.9%	46.1	46.8	4218	1241	0.07	0.42	246	266	57.3%	67.7%	39.1	45.6
ilEndoFla1	623	489	1.46	1.92	492	489	98.6%	99.2%	46.5	45.9	2110	1587	0.29	0.56	375	418	82.7%	89.0%	45.9	46.3
ilLymaMon1	251	301	10.31	5.77	917	912	99.2%	99.2%	46.5	46.8	5525	3162	0.17	0.47	610	633	64.6%	66.4%	41.6	43.8
ilNoctFim1	382	307	3.12	3.75	576	572	98.2%	98.9%	47.6	46.3	1493	1152	0.86	1.62	514	545	83.9%	94.9%	45.1	48.1
ilNotoUdd1	1128	809	1.52	2.07	829	814	98.9%	98.3%	45.5	46.8	4115	2791	0.26	0.69	629	780	75.3%	94.0%	42.8	45.6
ilParaStr1	218	129	4.42	6.92	480	481	99.4%	99.4%	48.0	47.4	1166	1309	0.64	1.78	413	431	84.3%	84.7%	47.1	46.9
ilRecuLeu1	1268	1247	1.09	0.98	748	746	98.5%	98.2%	44.4	45.1	4890	4029	0.12	0.23	416	534	54.1%	70.4%	42.0	43.4
ilThyaBat1	186	88	3.29	7.15	316	316	98.5%	98.9%	46.2	47.0	650	1073	0.89	2.77	296	327	91.6%	96.0%	46.9	45.3
ilVaneAta1	80	48	10.13	12.12	368	368	99.1%	99.3%	47.9	48.3	2185	1444	0.14	4.80	205	369	60.1%	97.7%	45.3	46.9
idAnopAqu88	260	79	18.22	15.05	188	181	99.2%	97.8%	49.5	50.6	1056	1833	0.11	4.02	90	211	43.2%	94.2%	48.9	43.8
idAnopCol22p13	272	189	5.48	5.01	260	260	99.4%	99.2%	49.4	49.5	1112	2474	0.12	0.23	77	148	27.1%	52.5%	43.9	40.3

	within 10%
	>10% better
	>25% better
	worse

RESULTS: BUG GENOMES

ALL SAMPLES: IPA vs. FALCON

Species	Primary										Haplotigs									
	# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV		# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA
ilAutoGam1	85	79	10.97	12.01	381	368	99.0%	99.2%	47.1	47.1	378	1613	6.98	9.09	330	379	93.8%	97.9%	46.3	44.9
ilCosmTra1	2173	1924	0.79	0.88	872	862	97.9%	97.5%	44.4	45.3	3932	3505	0.27	0.32	717	742	85.6%	89.8%	44.1	43.3
ilCranLig1	220	267	7.32	3.50	438	436	99.0%	98.9%	46.1	46.8	4218	1241	0.07	0.42	246	266	57.3%	67.7%	39.1	45.6
ilEndoFla1	623	489	1.46	1.92	492	489	98.6%	99.2%	46.5	45.9	2110	1587	0.29	0.56	375	418	82.7%	89.0%	45.9	46.3
ilLymaMon1	251	301	10.31	5.77	917	912	99.2%	99.2%	46.5	46.8	5525	3162	0.17	0.47	610	633	64.6%	66.4%	41.6	43.8
ilNoctFim1	382	307	3.12	3.75	576	572	98.2%	98.9%	47.6	46.3	1493	1152	0.86	1.62	514	545	83.9%	94.9%	45.1	48.1
ilNotoUdd1	1128	809	1.52	2.07	829	814	98.9%	98.3%	45.5	46.8	4115	2791	0.26	0.69	629	780	75.3%	94.0%	42.8	45.6
ilParaStr1	218	129	4.42	6.92	480	481	99.4%	99.4%	48.0	47.4	1166	1309	0.64	1.78	413	431	84.3%	84.7%	47.1	46.9
ilRecuLeu1	1268	1247	1.09	0.98	748	746	98.5%	98.2%	44.4	45.1	4890	4029	0.12	0.23	416	534	54.1%	70.4%	42.0	43.4
ilThyaBat1	186	88	3.29	7.15	316	316	98.5%	98.9%	46.2	47.0	650	1073	0.89	2.77	296	327	91.6%	96.0%	46.9	45.3
ilVaneAta1	80	48	10.13	12.12	368	368	99.1%	99.3%	47.9	48.3	2185	1444	0.14	4.80	205	369	60.1%	97.7%	45.3	46.9
idAnopAqu88	260	79	18.22	15.05	188	181	99.2%	97.8%	49.5	50.6	1056	1833	0.11	4.02	90	211	43.2%	94.2%	48.9	43.8
idAnopCol22p13	272	189	5.48	5.01	260	260	99.4%	99.2%	49.4	49.5	1112	2474	0.12	0.23	77	148	27.1%	52.5%	43.9	40.3

	within 10%
	>10% better
	>25% better
	worse

Very similar in size,
completeness & accuracy

RESULTS: BUG GENOMES

ALL SAMPLES: IPA vs. FALCON

Species	Primary # contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	Falcon	IPA	Falcon	IPA						
					Falcon	IPA	Falcon	IPA	Falcon	IPA
ilAutoGam1	85	79	10.97	12.01	381	368	99.0%	99.2%	47.1	47.1
ilCosmTra1	2173	1924	0.79	0.88	872	862	97.9%	97.5%	44.4	45.3
ilCranLig1	220	267	7.32	3.50	438	436	99.0%	98.9%	46.1	46.8
ilEndoFla1	623	489	1.46	1.92	492	489	98.6%	99.2%	46.5	45.9
ilLymaMon1	251	301	10.31	5.77	917	912	99.2%	99.2%	46.5	46.8
ilNoctFim1	382	307	3.12	3.75	576	572	98.2%	98.9%	47.6	46.3
ilNotoUdd1	1128	809	1.52	2.07	829	814	98.9%	98.3%	45.5	46.8
ilParaStr1	218	129	4.42	6.92	480	481	99.4%	99.4%	48.0	47.4
ilRecuLeu1	1268	1247	1.09	0.98	748	746	98.5%	98.2%	44.4	45.1
ilThyaBat1	186	88	3.29	7.15	316	316	98.5%	98.9%	46.2	47.0
ilVaneAta1	80	48	10.13	12.12	368	368	99.1%	99.3%	47.9	48.3
idAnopAqu88	260	79	18.22	15.05	188	181	99.2%	97.8%	49.5	50.6
idAnopCol22p13	272	189	5.48	5.01	260	260	99.4%	99.2%	49.4	49.5

Haplotigs									
# contigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA	Falcon	IPA
378	1613	6.98	9.09	330	379	93.8%	97.9%	46.3	44.9
3932	3505	0.27	0.32	717	742	85.6%	89.8%	44.1	43.3
4218	1241	0.07	0.42	246	266	57.3%	67.7%	39.1	45.6
2110	1587	0.29	0.56	375	418	82.7%	89.0%	45.9	46.3
5525	3162	0.17	0.47	610	633	64.6%	66.4%	41.6	43.8
1493	1152	0.86	1.62	514	545	83.9%	94.9%	45.1	48.1
4115	2791	0.26	0.69	629	780	75.3%	94.0%	42.8	45.6
1166	1309	0.64	1.78	413	431	84.3%	84.7%	47.1	46.9
4890	4029	0.12	0.23	416	534	54.1%	70.4%	42.0	43.4
650	1073	0.89	2.77	296	327	91.6%	96.0%	46.9	45.3
2185	1444	0.14	4.80	205	369	60.1%	97.7%	45.3	46.9
1056	1833	0.11	4.02	90	211	43.2%	94.2%	48.9	43.8
1112	2474	0.12	0.23	77	148	27.1%	52.5%	43.9	40.3

	within 10%
	>10% better
	>25% better
	worse

Very similar in size,
completeness & accuracy


	Improved haplotype separation
	>90% haplotype resolved

RESULTS: BUG GENOMES


ALL SAMPLES: IPA vs. HICANU¹

Species	Primary		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	# contigs									
	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA
ilAutoGam1	128	79	12.17	12.01	375	368	99.1%	99.2%	47.1	47.1
ilCosmTra1		1924		0.88		862		97.5%		45.3
ilCranLig1	211	267	4.94	3.50	435	436	99.1%	98.9%	48.7	46.8
ilEndoFla1	662	489	1.56	1.92	488	489	99.3%	99.2%	46.8	45.9
ilLymaMon1	157	301	13.59	5.77	912	912	99.2%	99.2%	47.3	46.8
ilNoctFim1	783	307	1.88	3.75	577	572	98.8%	98.9%	48.2	46.3
ilNotoUdd1	1147	809	1.56	2.07	826	814	99.1%	98.3%	46.1	46.8
ilParaStr1	105	129	11.77	6.92	482	481	99.1%	99.4%	48.3	47.4
ilRecuLeu1		1247		0.98		746		98.2%		45.1
ilThyaBat1	214	88	3.31	7.15	319	316	98.4%	98.9%	44.3	47.0
ilVaneAta1	242	48	12.18	12.12	372	368	99.1%	99.3%	48.6	48.3
idAnopAqu88		79		15.05		181		97.8%		50.6
idAnopCol22p13		189		5.01		260		99.2%		49.5

Species	Haplotigs		Contig N50 (Mb)		Size (Mb)		BUSCO C		QV	
	# contigs									
	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA	HiCanu	IPA
ilAutoGam1	643	1613	8.14	9.09	356	379	97.0%	97.9%	46.5	44.9
ilCosmTra1		3505		0.32		742		89.8%		43.3
ilCranLig1	2268	1241	0.36	0.42	425	266	95.2%	67.7%	46.6	45.6
ilEndoFla1	2155	1587	0.39	0.56	442	418	95.0%	89.0%	46.1	46.3
ilLymaMon1	4208	3162	0.57	0.47	950	633	96.3%	66.4%	45.7	43.8
ilNoctFim1	1258	1152	0.69	1.62	530	545	96.8%	94.9%	47.6	48.1
ilNotoUdd1	3193	2791	0.49	0.69	802	780	95.4%	94.0%	46.2	45.6
ilParaStr1	789	1309	1.94	1.78	465	431	95.5%	84.7%	47.5	46.9
ilRecuLeu1		4029		0.23		534		70.4%		43.4
ilThyaBat1	900	1073	0.81	2.77	318	327	97.4%	96.0%	43.7	45.3
ilVaneAta1	769	1444	4.31	4.80	357	369	97.7%	97.7%	47.8	46.9
idAnopAqu88		1833		4.02		211		94.2%		43.8
idAnopCol22p13		2474		0.23		148		52.5%		40.3



Very similar in size,
completeness & accuracy




¹<https://github.com/darwintreeoflife/darwintreeoflife.data/>

UPCOMING FEATURES

UPCOMING FEATURES

—Integration of “purge_dups” into the workflow



—Phasing improvements

—Read tracking improvements for better polishing

AVAILABILITY, INSTALLATION AND USAGE

AVAILABILITY

- IPA available on Bioconda!
- More details and documentation available here:
 - <https://github.com/PacificBiosciences/pbbioconda/wiki/Improved-Phased-Assembler>
 - <https://github.com/PacificBiosciences/pbipa>
 - <https://github.com/PacificBiosciences/pbbioconda>



Ivan Sovic @IvanSovic · May 28


Proud to announce the team @PacBio and myself are working on a new Improved and Phased Assembly method for HiFi reads called IPA!

Fast, contiguous, runs locally and on a cluster!

Early version now on Bioconda, package "pbipa".

[github.com/PacificBioscie...](https://github.com/PacificBiosciences/pbbioconda/wiki/Improved-Phased-Assembler)

@zevkronenberg @drsarahdoom



6

66

130



INSTALLATION AND USAGE

—Installation

```
conda create -n ipa -c bioconda -c conda-forge -c defaults
conda activate ipa
conda install pbipa
```

—Run assembly on a local machine:

```
ipa local --nthreads 24 --njobs 1 -i <reads.fasta>
```

—Run assembly on an SGE cluster:

```
ipa dist --nthreads 24 --njobs 40 -i <reads.fasta> \
    --cluster-args 'qsub -S /bin/bash -N ipa.{rule} -cwd -q default -pe smp {params.num_threads} -e
qsub_log/ -o qsub_log/ -V'
```

SUMMARY

- IPA delivers highly accurate and contiguous assemblies, with high speed and accurately phased haplotig regions
 - Generates true haplotigs constructed through a phasing process
- Polishes the phased genome to achieve **>Q50** accuracy!
- Further evaluations and developments ongoing
- Potential for IPA & HiCanu to learn from each other
- Ease of use!
- **Work in progress:**
 - Integrate "purge_dups" directly into the workflow
 - Improve contiguity of the phased assembly
 - Optimization of all stages

THANK YOU!



IPA TEAM


Ivan Sović
Zev Kronenberg
Christopher Dunn
Sarah Kingan
Derek Barnett
James Drake
Jonas Korf

PACBIO

Armin Töpfer
Paul Peluso
Greg Conception

COLLABORATORS


Jay Ghurye
Nathan Truelove
Barbara Block
Mara Lawniczak
Darwin Tree of Life Project



PUBLIC HIFI DATA

HG002 Human Pan-Genome Reference Consortium

- 4 cells: 2 cells 20kb and 2 cells 15kbp
- ~34x coverage
- https://github.com/human-pangenomics/HG002_Data_Freeze_v1.0



human-pangenomics / HG002_Data_Freeze_v1.0

Watch 6 Star 26 Fork 3

Code Issues 0 Pull requests 0 Actions Projects 0 Wiki Security Insights

Human Pangenome Reference Consortium - HG002 Data Freeze (v1.0)

32 commits 1 branch 0 packages 0 releases 2 contributors

Branch: master New pull request Create new file Upload files Find file Clone or download

Sequencing Data

The annotated table of sequence data can be downloaded [here](#).

HG002 Data Freeze (v1.0) Recommended downsampled data mix

We encourage assembly groups to use as much of the data from the HG002 freeze as possible to get the best assembly they can. However, as no two groups are likely to use exactly the same subset of data, making comparison more difficult, and the size and variety of the HG002 freeze is not representative of what is likely to be available in future freezes, we recommend that assembly groups also run their pipeline on the following set of 4 downsampled datasets from the HG002 (NA24385) human cell line:

PacBio HiFi:

~34X coverage of Sequel II System with Chemistry 2.0

15kb:

- https://s3-us-west-2.amazonaws.com/human-pangenomics/HG002/hpp_HG002_NA24385_son_v1/PacBio_HiFi/15kb/m64012_190920_173625.Q20.fastq
- https://s3-us-west-2.amazonaws.com/human-pangenomics/HG002/hpp_HG002_NA24385_son_v1/PacBio_HiFi/15kb/m64012_190921_234837.Q20.fastq

20kb:

- https://s3-us-west-2.amazonaws.com/human-pangenomics/HG002/hpp_HG002_NA24385_son_v1/PacBio_HiFi/20kb/m64011_190830_220126.Q20.fastq
- https://s3-us-west-2.amazonaws.com/human-pangenomics/HG002/hpp_HG002_NA24385_son_v1/PacBio_HiFi/20kb/m64011_190901_095311.Q20.fastq

PUBLIC HIFI DATA

CHM13 data from the HiCanu preprint

— 5 HiFi datasets

— <https://www.ncbi.nlm.nih.gov/sra/?term=PRJNA530776>

☐ [WGS of CHM13 with PacBio CCS](#)

1. 1 PACBIO_SMRT (Sequel II) run: 1M spots, 21G bases, 15.7Gb downloads
Accession: SRX7897688

☐ [WGS of CHM13 with PacBio CCS](#)

2. 1 PACBIO_SMRT (Sequel II) run: 1.4M spots, 28.7G bases, 21.7Gb downloads
Accession: SRX7897687

☐ [WGS of CHM13 with PacBio CCS](#)

3. 1 PACBIO_SMRT (Sequel II) run: 1.6M spots, 25.6G bases, 16.3Gb downloads
Accession: SRX7897686

☐ [WGS of CHM13 with PacBio CCS](#)

4. 1 PACBIO_SMRT (Sequel II) run: 1.6M spots, 25.1G bases, 16Gb downloads
Accession: SRX7897685

☐ [WGS of CHM13 with PacBio CCS](#)

5. 4 PACBIO_SMRT (Sequel II) runs: 6.9M spots, 75.6G bases, 47.3Gb downloads
Accession: SRX5633451



The screenshot shows the bioRxiv preprint interface. At the top, there's a navigation bar with links: HOME | ABOUT | SUBMIT | NEWS & NOTES | ALERTS / RSS | CHANNELS. Below this is a search bar with the text "Search" and a magnifying glass icon. To the right of the search bar is a link for "Advanced Search".

The main content area features a yellow banner with the text: "bioRxiv is receiving many new papers on coronavirus SARS-CoV-2. A reminder: these are preliminary reports that have not been peer-reviewed. They should not be regarded as conclusive, guide clinical practice/health-related behavior, or be reported in news media as established information."


Below the banner, the preprint title "HiCanu: accurate assembly of segmental duplications, satellites, and allelic variants from high-fidelity long reads" is displayed. The authors listed are: Sergey Nurk, Brian P. Walenz, Arang Rhie, Mitchell R. Vollger, Glennis A. Logsdon, Robert Grothe, Karen H. Miga, Evan E. Eichler, Adam M. Phillippy, and Sergey Koren. The DOI is https://doi.org/10.1101/2020.03.14.992248. A note states: "This article is a preprint and has not been certified by peer review [what does this mean?]."

On the left side of the preprint, there are links for "New Results", "Comment on this paper", "Previous", and "Next".

On the right side, there are links for "Download PDF", "Supplementary Material", "Data/Code", "XML", and "Revision Summary". There are also social media links for "Email", "Share", and "Citation Tools".

At the bottom, there's a navigation bar with links: Abstract, Full Text, Info/History, Metrics, and a "Preview PDF" button.

PUBLIC HIFI DATA



HG002

15 kb + 20 kb library

6 SMRT Cell 8M

[Data: PRJNA586863](https://www.ncbi.nlm.nih.gov/bioproject/PRJNA586863)




Oryza sativa indica MH63

17 kb + 24 kb library

2 SMRT Cell 8M

[Data: PRJNA573706](https://www.ncbi.nlm.nih.gov/bioproject/PRJNA573706)



Drosophila melanogaster F1

19 kb + 24 kb library

2 SMRT Cell 8M

[Data: PRJNA573706](https://www.ncbi.nlm.nih.gov/bioproject/PRJNA573706)










PUBLIC HIFI DATA



bioRxiv
THE PREPRINT SERVER FOR BIOLOGY

New Results

Highly accurate long-read HiFi sequencing data for five complex genomes

Ting Hon, Kristin Mars, Greg Young,  Yu-Chih Tsai, Joseph W. Karalius, Jane M. Landolin,  Nicholas Maurer,
 David Kudrna, Michael A. Hardigan,  Cynthia C. Steiner,  Steven J. Knapp,  Doreen Ware,
 Beth Shapiro,  Paul Peluso,  David R Rank

doi: <https://doi.org/10.1101/2020.05.04.077180>





www.pacb.com

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.