



PACIFIC  
BIOSCIENCES®




# SMRT Sequencing – How it Works

Elizabeth Tseng, Principal Scientist, PacBio


 @Magdoll

# SEQUENCE WITH CONFIDENCE


PacBio long-read sequencing provides access to the full spectrum of genetic variation, driving discovery across all fields of life science.




PacBio Systems are powered by **Single Molecule, Real-Time (SMRT) Sequencing**, a technology proven to produce exceptionally long reads with high accuracy.



# FROM SAMPLE TO SMRT SEQUENCING




# OUR CORE TECHNOLOGY






# GENERATE HIGHLY ACCURATE LONG READS




# PACBIO SEQUENCING 101 CLIP




[PacBio: Sequencing 101](#)

# SMRT SEQUENCING ENABLES THE FULL RANGE OF PACBIO APPLICATIONS



## Accelerate Your Science







# The Iso-Seq Method

From Sample Prep to Analysis

# TRADITIONAL RNA-SEQ



# ISO-SEQ METHOD = FULL-LENGTH TRANSCRIPT SEQUENCING




No assembly required



## Iso-Seq Express kit

- 60-300 ng total RNA
- Full-length cDNA
- Multiplexing support





	FL Reads	Unique Genes	Unique Transcripts
UHRR	4,734,362	16,328	183,689
Alzheimer Brain	4,277,293	17,670	162,290

[Dataset: UHRR](#)

[Dataset: Alzheimer Brain](#)



SMRT  
SEQUENCING


**1 DAY**




## Sequel II System

- 1 SMRT Cell 8M for whole transcriptome
- Up to 4 million full-length reads

# FULL BIOINFORMATICS SOLUTION




# CLASSIFICATION OF TRANSCRIPTS



# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- **Comprehensive**
- Full-Length
- Highly accurate



162,290 transcripts

80 – 14,288 bp  
(mean: 3.3 kb)

[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- **Comprehensive**
- Full-Length
- Highly accurate

	Known	Novel	Total
<b>Genes</b>	17,051	619	<b>17,670</b>
<b>Isoforms</b>	51,660	110,630	<b>162,290</b>

[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- Comprehensive
- **Full-Length**
- Highly accurate

Category	Count	Description
<b>FSM</b>	32,649	Perfect match
<b>ISM</b>	19,011	Incomplete match
<b>NIC</b>	84,610	Novel isoform using known junctions
<b>NNC</b>	25,323	Novel isoform using at least novel junction
<b>Antisense</b>	321	Anti-sense to known gene
<b>Intergenic</b>	376	Intergenic

[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- Comprehensive
- **Full-Length**
- Highly accurate

Category	Count	CAGE peak within 50 bp
<b>FSM</b>	32,649	70%
<b>ISM</b>	19,011	37%
<b>NIC</b>	84,610	36%
<b>NNC</b>	25,323	57%
<b>Antisense</b>	321	24%
<b>Intergenic</b>	376	24%

[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- Comprehensive
- **Full-Length**
- Highly accurate


Category	Count	CAGE peak within 50 bp	polyA Motif Detected
<b>FSM</b>	32,649	70%	72%
<b>ISM</b>	19,011	37%	62%
<b>NIC</b>	84,610	36%	55%
<b>NNC</b>	25,323	57%	72%
<b>Antisense</b>	321	24%	43%
<b>Intergenic</b>	376	24%	38%

[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)



# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM


- Comprehensive
- **Full-Length**
- Highly accurate



[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM


- Comprehensive
- **Full-Length**
- Highly accurate



[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM


- Comprehensive
- Full-Length
- **Highly accurate**



[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

# ISO-SEQ METHOD ON THE SEQUEL II SYSTEM

- Comprehensive
- Full-Length
- **Highly accurate**



[Dataset: Alzheimer brain on 1 SMRT Cell 8M](#)

## Iso-Seq Express kit

- 60-300 ng total RNA
- Full-length cDNA
- Multiplexing support



LIBRARY  
PREP

**1 DAY**



SMRT  
SEQUENCING

**1 DAY**



DATA  
ANALYSIS

**1 DAY**

## Bioinformatics

- Iso-Seq3 in SMRT Analysis
- Reads to ORF in 1 day
- Downstream community tools



## Sequel II System

- 1 SMRT Cell 8M for whole transcriptome
- 4 million full-length reads



# Research Highlights

# GENETIC DIAGNOSIS FOR RARE DISEASES

Circulation: Genomic and Precision Medicine

## RESEARCH LETTER


### Targeted Long-Read RNA Sequencing Demonstrates Transcriptional Diversity Driven by Splice-Site Variation in *MYBPC3*

To date, clinical sequencing has focused on genomic DNA using targeted panels and exome sequencing. Sequencing of a large hypertrophic cardiomyopathy (HCM) cohort revealed that positive identification of a disease-associated variant was returned in only 32% of patients, with an additional 15% receiving inconclusive results.<sup>1</sup> When genome sequencing fails to reveal causative variants, the transcriptome may provide additional diagnostic clarity. A recent study examining patients with genetically undiagnosed muscle disorders found that RNA sequencing, when used as a complement to exome and whole genome sequencing, had an overall diagnosis rate of 35%.<sup>2</sup>

Alexandra Dainis, PhD  
Elizabeth Tseng, PhD  
Tyson A. Clark, PhD  
Ting Hon, MS  
Matthew Wheeler, MD, PhD  
Euan Ashley, MB, ChB, DPhil

<https://medium.com/@magdoll>

ing in PacBio  
7 min read



### Finding Diagnosis in a Sea of Transcripts: The Case of a Hypertrophic Cardiomyopathy Patient

# GENETIC DIAGNOSIS FOR RARE DISEASES

- 21y female patient with severe hypertrophic cardiomyopathy (HCM)
- HCM panel identified single base (c.1898-1G>A) mutation in MYBPC3
- Mutation expected to affect splicing between exon 19-20

Circulation: Genomic and Precision Medicine

## RESEARCH LETTER



### Targeted Long-Read RNA Sequencing Demonstrates Transcriptional Diversity Driven by Splice-Site Variation in *MYBPC3*

**T**o date, clinical sequencing has focused on genomic DNA using targeted panels and exome sequencing. Sequencing of a large hypertrophic cardiomyopathy (HCM) cohort revealed that positive identification of a disease-associated variant was returned in only 32% of patients, with an additional 15% receiving inconclusive results.<sup>1</sup> When genome sequencing fails to reveal causative variants, the transcriptome may provide additional diagnostic clarity. A recent study examining patients with genetically undiagnosed muscle disorders found that RNA sequencing, when used as a complement to exome and whole genome sequencing, had an overall diagnosis rate of 35%.<sup>2</sup>



Alexandra Dainis, PhD  
Elizabeth Tseng, PhD  
Tyson A. Clark, PhD  
Ting Hon, MS  
Matthew Wheeler, MD, PhD  
Euan Ashley, MB, ChB, DPhil



# GENETIC DIAGNOSIS FOR RARE DISEASES





# GENETIC DIAGNOSIS FOR RARE DISEASES



Dainis et al., “Targeted Long-Read RNA Sequencing Demonstrates Transcriptional Diversity Driven by Splice-Site Variation in MYBPC3.” *Circulation. Genomic and Precision Medicine* (2019)

# GENETIC DIAGNOSIS FOR RARE DISEASES




Dainis et al., “Targeted Long-Read RNA Sequencing Demonstrates Transcriptional Diversity Driven by Splice-Site Variation in MYBPC3.” *Circulation. Genomic and Precision Medicine* (2019)


## GENETIC DIAGNOSIS FOR RARE DISEASES

- 21y female patient with severe hypertrophic cardiomyopathy (HCM)
- HCM panel identified single base (c.1898-1G>A) mutation in MYBPC3
- Mutation expected to affect splicing between exon 19-20
- Iso-Seq method identified patient-specific alternative splice variants (AS 1.1-1.11) linked to the c.1898-1G>A mutation that are not expressed on the WT allele

# FUSION TRANSCRIPTS IN CANCER




# FUSION TRANSCRIPTS IN CANCER




- In B-cell development, one IGH allele is silenced
- IGH-DUX4 translocation is subtype of B-cell acute lymphoblastic leukemia (B-ALL)
- Accounts for ~7% of pediatric B-ALL
- DUX4 is a transcription factor located in a GC-rich region
- Short reads difficult to map accurately IGH-DUX4 translocation

# ISO-SEQ SEQUENCING OF NALM6 B-ALL CELL LINE




# FUSION TRANSCRIPTS IN CANCER




- IGHM locus is diploid with three SNPs (A-A-A vs T-G-G)
- Iso-Seq analysis shows functional IgM comes exclusively from A-A-A wildtype
- Iso-Seq analysis shows antisense DUX4 comes exclusively from T-G-G fusion allele --> proves that IGH-DUX4 translation comes from the silenced IGH allele





# ISO-SEQ METHOD FOR PLANT & ANIMAL



## Brahman x Angus F1 Cattle


- Allele-specific isoform expression
- Tissue-specific isoform expression

[Low et al \(2020\)](#)

## Cannabis

- Tissue-specific transcripts associated w/ THC & CBD synthesis
- Chr Y gene annotation

McKernan, [webinar](#)



## Grizzly Bear


- Tissue-specific alternative splicing
- Hibernation vs active state

Trojahn, Kelley, et al. (WSU)

# LEARN MORE ABOUT AND THE ISO-SEQ METHOD



[pacb.com/applications/rna-sequencing/](https://pacb.com/applications/rna-sequencing/)



[tinyurl.com/isoseq-google](https://tinyurl.com/isoseq-google)



@PacBio @Magdoll



[github.com/PacificBiosciences](https://github.com/PacificBiosciences)  
[github.com/Magdoll](https://github.com/Magdoll)



PACBIO<sup>®</sup>

[www.pacb.com](http://www.pacb.com)

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2020 by Pacific Biosciences of California, Inc. All rights reserved. Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT, SMRTbell, Iso-Seq, and Sequel are trademarks of Pacific Biosciences. Pacific Biosciences does not sell a kit for carrying out the overall No-Amp Targeted Sequencing method. Use of these No-Amp methods may require rights to third-party owned intellectual property. FEMTO Pulse and Fragment Analyzer are trademarks of Agilent Technologies Inc.

All other trademarks are the sole property of their respective owners.