


Single Cell Sample Preparation



Diana Burkart-Waco, PhD

dburkart@ucdavis.edu

SRA at UC Davis DNATECH Core


<https://dnatech.genomecenter.ucdavis.edu/>




Why single cell?

- Methods – SC isolation
 - Single cell isolation
 - Single nuclei
- Methods - sample QC
 - Do I have single cells?
 - Are they alive?
 - Are they too big?
 - Did I isolate the correct cells?
- Cell labeling
- Single cell studies at DNA Tech


scRNA-seq



10X Genomics



- Captures 100-10,000+ cells per sample.
- Recovers up to ~65% of cells.
- Low doublet rate.
– (~0.9% per 1,000 cells).
- 40 micron size limit.
- Expensive!




Alternative technologies

- Homemade drop-seq.
- Qiagen UPX.
- BD Rhapsody.
- SPLIT-seq.
- Dolomite Nadia System.
- Fluidigm.
- And the list goes on...


Need both **single cells** and to maintain **sample composition**.



SPLIT-Seq



- Use combinatorial barcoding to label cells of origin.
- Great because don't need fancy equipment.
- Increased cell size limits.
- Input = fixed cells or nuclei.

<https://www.splitbiosciences.com/>





- Inject cells into cartridge and fall into microwells with barcoded beads.
- Similar cell size limit and cost to 10X.
- Higher capture rate.
- More flexibility in targeted capture with custom beads.

- Why single cell?



Methods – SC isolation

- Single cell isolation
 - Single nuclei
- Methods - sample QC
 - Do I have single cells?
 - Are they alive?
 - Are they too big?
 - Did I isolate the correct cells?
 - Cell labeling
 - Single cell studies at DNA Tech

- Cell isolation guides available at:

<https://www.support.10xgenomics.com/single-cell-gene-expression/sample-prep/>


Sample Prep

▼ Demonstrated Protocol (14 documents)

- [Single Cell Gene Expression Demonstrated Protocol Compatibility Table](#)
- [Cell Surface Protein Labeling for Single Cell RNA Sequencing Protocols](#)
- [Methanol Fixation of Cells for Single Cell RNA Sequencing](#)
- [Isolation of Nuclei for Single Cell RNA Sequencing](#)
-  [Single Cell Protocols - Cell Preparation Guide](#)
- [Enrichment of CD3+ T Cells from Dissociated Tissues for Single Cell RNA Sequencing and Immune Repertoire Profiling](#)
- [Tumor Dissociation for Single Cell RNA Sequencing](#)
- [Thawing Dissociated Tumor Cells for Single Cell RNA Sequencing](#)
- [Single Cell Suspensions from Cultured Cell Lines for Single Cell RNA Sequencing](#)
- [Removal of Dead Cells from Single Cell Suspensions for Single Cell RNA Sequencing](#)
- [Moss Protoplast Suspension for Single Cell RNA Sequencing](#)
- [Fresh Frozen Human-Mouse Cell Line Mixtures for Single Cell RNA Sequencing](#)
- [Fresh Frozen Human Peripheral Blood Mononuclear Cells for Single Cell RNA Sequencing](#)
- [Dissociation of Mouse Embryonic Neural Tissue for Single Cell RNA Sequencing](#)



- Cell isolation → biggest source technical variation.
- Dissociation and preparation depends on cell type.

General Workflow



Tissue Dissociation


- Many ways to dissociate tissue:
 - Mechanical
 - Enzymatic
 - Automated blending
 - Microfluidics devices
- Considerations:
 - Speed
 - Consistency in results
 - Good representation of all cell types



Arabidopsis
protoplasts

Same method,
different days.

Factors Influencing Dissociation



Unfortunately, SC sample optimization is best achieved through trial and error...


Single cell vs. nuclei

- Single cell captures more transcripts, but is a harder protocol.
- When to use nuclei:
 - Cells cannot be harvested intact or viable (e.g. adipocytes, neurons).
 - Cells are too big for capture (e.g. cardiomyocytes).
 - Tissue frozen.

Tissue dissection



Cell lysis



Detergent
Grinding

Nuclei purification



Flow cytometry
Density Gradient

Many protocols available...

- Start here:

Customer Developed Protocol

'Frankenstein' protocol for nuclei isolation from fresh and frozen tissue



Contributed by:
Luciano Martelotto, Ph.D., Melbourne, Centre for Cancer
Comprehensive Cancer Centre




www.collaslab.com

Isolation of Nuclei from Somatic Cells

1. HeLa Cells, 293T Cells, NT2 Cells

Cell preparation

- harvest cells from flasks as per standard protocol
- spin cells in 50 ml conical tube at 1500 rpm for 10 min at RT
- resuspend cells in 30 ml PBS; take a 50 μ l sample to determine concentration



10x Genomics®

Sample Preparation Demonstrated Protocols

Isolation of Nuclei for Single Cell RNA
Sequencing

- Why single cell?
- Methods – SC isolation
 - Single cell isolation
 - Single nuclei



Methods - sample QC


- Do I have single cells?
- Are they alive?
- Are they too big?
- Did I isolate the correct cells?
- Cell labeling
- Single cell studies at DNA Tech

Factors influencing success


- Cell debris / dead cells
 - Clog microfluidics, free floating RNA → noise.
- Aggregates
 - No longer single cell data.
 - Clog microfluidics.
- Buffer
 - Inhibit polymerase → decrease library complexity.
- Storage conditions

But the most important factor is **cell / nuclei counting!**

Workflow



Tissue aggregates



Tissue clump



Single nucleus

DAPI stain


60X

No scale bar

Recommended treatment: optimize tissue dissociation

Cell debris

- Mouse DRG.



Clean DRG sample



Noisy sample

Recommended treatment: filtration, centrifuge, add blocking agent.

- Lettuce nuclei prep.






Photo credit: Mohan Prem Anand Marimuthu

Recommended treatment: modify dissociation times / detergent concentrations, change density gradient.

Doublets

- Non-single cell clumps.
- Integrated into droplets and cannot* distinguish from single cells.




Bad for any
single cell
experiment.



Recommended treatment: filtration, change [blocking agent]

Dead cells (debris?)



Recommended treatment: dead cell removal.


<https://www.miltenyibiotec.com/US-en/products/dead-cell-removal-kit.html>

Nuclei QC


- Automated cell counters suboptimal, but possible with fluorescence.
 - Countess II FL, Devonix CellDrop.
- Stains
 - Trypan blue → 100% dead (not great).
 - Ethidium homodimer stain (good for excluding debris).
 - DAPI, PI other options




60X: DAPI



Hemocytometer




CountessII




Credit: 10X

Workflow



Automated counting

- Countess II.



- Pros (+):



- Fast.
- Live/dead cell counts.
- Cell size estimates.

- Cons (-):

- Cell size limits (4-30 μm).
- Doesn't do well with odd cell shapes, debris, nuclei.


Manual counting

- Hemocytometer.



- Pros (+):
 - Reliable cell counts.
 - Count small cells.
 - Visualize cells vs. debris.
- Cons (-):
 - Slow...





- FACS






<https://flowcytometry-embl.de/cell-sorting/>

- Pros (+):
 - Sort live / dead.
 - Enrich cells of interest.
 - Determine whether correct cells isolated (qPCR works too).
- Cons (-):
 - Overestimates cell counts (~2x).

The good...




The bad...



Challenges Cell Recovery

- 10X is optimized for model cell types – human PBMCs, etc.
 - Round, easy to count, size well below size of microfluidics.
 - These yield consistent results. 10K target. 8K-10K recovered.



Estimated Number of Cells

8,266


Mean Reads per Cell

21,128

Median Genes per Cell

1,027

- But most experiments outside of culture don't look like this.



Variable size,
shape, and
viability

How maximize output?

- Carefully craft **experimental design** and **sample prep** .
- QC **cells** before real sample set-up.
- Concentration: aim for the median.
 - **700-1,200 cells per μ l optimal.**
 - Too high \rightarrow dilute.
 - Too low \rightarrow tough one (concentration impacts yield, pooling replicates suboptimal).
 - Count in replicates (at least $n=2$).
- Viability: **70% minimum.**
 - Nuclei and methanol fixed cells (0%).
- Treat cells gently.
 - Wide bore pipette tips, keep cells at preferred temp, work quickly..

- Why single cell?
- Methods – SC isolation
 - Single cell isolation
 - Single nuclei
- Methods - sample QC
 - Do I have single cells?
 - Are they alive?
 - Are they too big?
 - Did I isolate the correct cells?





Cell labeling

- Single cell studies at DNA Tech


Cell Hashing

- Label cells (either oligo-tagged antibodies) or oligos.
- Allows for multiplexing and super-loading scRNA-seq platforms – cost saving and minimize tech variation.
- Many options available:
 - MULTI-seq (C. McGinnis / Z Gartner @ UCSF)
 - Biolegend TotalSeq
 - BD Single-Cell Multiplexing Kit
 - 10X CellPlex (coming later this year)




MultiSeq Tag Design

- Oligo targeted to plasma membrane with 5' lignoceric acid amide.
 - Co-anchor increases stability and specificity.
- 8bp barcode conjugated to 3' Poly-A capture sequence.



- Once in droplets, PolyA mimics endogenous transcripts and hybridizes to dT oligos in beads.



Barcode --> tag suspension

AAAs → capture transcripts

“MultiSeq” (C. McGinnis)


- Why single cell?
- Methods – SC isolation
 - Single cell isolation
 - Single nuclei
- Methods - sample QC
 - Do I have single cells?
 - Are they alive?
 - Are they too big?
 - Did I isolate the correct cells?
- Cell labeling



Single cell studies at DNA Tech

Assay summary (organism)

- 78 projects since July 2017 (primarily 3' GEX).



Sample distribution by organism



SCIENTIFIC REPORTS



OPEN

Molecular profiling of resident and infiltrating mononuclear phagocytes during rapid adult retinal degeneration using single-cell RNA sequencing

Kaitryn E. Ronning¹, Sarah J. Karlen², Eric B. Miller¹ & Marie E. Burns^{1,2,3}

Received: 15 October 2018

Accepted: 27 February 2019

Published online: 19 March 2019

CellPress

Sneak Peek

A PREVIEW OF PAPERS UNDER REVIEW

Single-Cell RNA-Seq Reveals Profound Alterations in Mechanosensitive Dorsal Root Ganglion Neurons with Vitamin E Deficiency

iScience

50 Pages • Posted: 5 Aug 2019 • Sneak Peek Status: **Review Complete**



[Carrie Finno](#)

University of California, Davis, School of Veterinary Medicine, Department of Population Health and Reproduction

[Janel Peterson](#)

University of California, Davis, School of Veterinary Medicine, Department of Population Health and Reproduction

'564' – challenges and solutions



- DRG fragile, difficult to count, odd sizes (clog 10X chip?).
- Count manually → better counts.
- Accepted some degree of induced technical variation
 - Prepped samples over three days.
 - Improved sample quality and reproducibility, introduced prep day variation.
 - As Matt mentioned mitigate with experimental design.

Important resources

- 10X Genomics
 - <https://support.10xgenomics.com/single-cell-gene-expression>
 - Paul Scott, Sales Executive (paul.scott@10xgenomics.com)
- UC Davis Flow Cytometry
 - http://www.ucdmc.ucdavis.edu/pathology/research/research_labs/flow_cytometry/index.html
 - Bridget McLaughlin (Technical Director)
- UC Davis DNA Technology Core
 - <http://dnatech.genomecenter.ucdavis.edu/single-cell-analyses/>
- UC Davis Bioinformatics Core
 - <https://bioinformatics.ucdavis.edu/>



Single Cell Expression Profiling & Genomics (10X Genomics, Fluidigm & Plate-based scRNA-Seq)

10X Genomics Single Cell Sequencing

Please schedule any 10X single cell experiment at least a week in advance. We highly advise a consultation prior to experiment scheduling.

Update March 2019: MULTI-Seq reagents promise to enable the labeling and pooling of up to 96 cell suspension samples. Reagents are in stock. [Please inquire!](#)

The 10X Genomics Single Cell system is the single-cell expression profiling platform enabling the analysis of large cell numbers at the highest capture efficiency (of up to 65%). The technology allows for high-throughput single cell transcriptomics of many different cell types as well as **single-nuclei expression profiling**. The flexible workflow encapsulates 500 to 20,000 cells or nuclei per library together with micro-beads into nano-droplets. Each bead is loaded with adapters containing one of 750,000 different barcodes for the single cell RNA-seq library preps. In contrast to other protocols (e.g. Drop-Seq) the 10X controller is capable of loading "all" droplets with micro-beads, enabling single-Poisson distribution loading and thus high capture efficiencies (in contrast to double-Poisson loading of other protocols). The single-cell encapsulating process is significantly faster compared to microfluidic or Drop-Seq. Up to eight samples can be processed per batch within minutes. The resulting data can be analyzed with the free Cell Ranger and Loupe Cell Browser software. In addition the Bioinformatics Core has developed a custom single-cell data analysis pipeline for 10X data.



10X Chromium Single Cell Features:

- Fast workflow from cell suspension to 3'-DNA library
- Captures 100-10,000+ cells per sample and up to 80,000 cells per run in < 7 minutes.
- Recovers up to ~65% of cells (typically 50%).

Cost 10X Reagents

Product Code	Product	List Price	Discount
1000120	Chromium Next GEM Chip G Single Cell Kit, 48 rxns	USD 1,480.00	
1000121	Chromium Next GEM Single Cell 3' GEM, Library & Gel Bead Kit v3.1, 16 rxns	USD 21,600.00	
		Quote Created Date	3/10/2020

- Chips are single use (n=1 to n=8 samples) → \$247 EACH
- Reagents (n=1) → \$1,350 EACH.
- Total (n=1, reagents only): \$1,597.
 - This doesn't include you or your cores labor!
- Sequencing? Bioinformatics?
 - Depends on many factors.
 - We charge \$1,000 for 350M reads (~10K cells with 35K reads per cell).

Acknowledgements

- UC Davis DNA Technology Core team.
- UC Davis Bioinformatics Core team.
- 10X Genomics.
 - Adam Bemis (Field Applications Scientist)
 - Paul Scott (Sales Executive)

